



VIDAS[®] ToRC PANEL

IgG - IgM - AVIDITY
3 markers for optimal safety



PIONEERING DIAGNOSTICS

Did you know?

- IgG is widely used to determine the presence of specific immunity.

Unlike false negative results, false positives or IgG over-sensitivity can have severe consequences.

- If IgM appear earlier than IgG, they also disappear earlier. This makes IgM an early and sensitive marker of acute infection.

The first quality of an IgM test is sensitivity, but IgM can persist for over 6 months and sometimes even for years.

- Serology is not always clear-cut. Reactivation, reinfection, polyclonal stimulation, residual IgM, non-specific IgM can lead to unnecessary treatment, amniocentesis, abortion and distress.^(1,2)

The Solution = Avidity

Right from the first sample collection, avidity measurement excludes most recent infections.



Serodiagnosis of ToRC in pregnant women

3 infections


Toxoplasmosis, Rubella and Cytomegalovirus can lead to severe sequelae for future newborns.

	Fq of C.I.*
CMV	1 newborn / 100 to 1,000
Toxo	1 newborn / 1,000 to 10,000
Rub	1 newborn / 10,000 to 100,000

* Frequency of congenital infection industrialized countries^(3,5)

3 markers

Diagnosis of ToRC is based on serology and comprises 3 markers: IgG, IgM, IgG Avidity.



3 test steps

Step 1	Step 2	Step 3
Immunity status	Suspected infection	Exclusion of recent infection
IgG detection	IgM detection	Avidity measurement

VIDAS® ToRC PANEL

IgG - IgM - AVIDITY 3 markers for optimal safety

With VIDAS® solution:

- Easy-to-use platform
- Choose mother and child safety
- The 3 markers answer all diagnostic guidelines for pregnancy follow-up

A

Immunity status ToRC IgG

Case study: Evaluation of a Toxo IgG assay⁽⁶⁾

sample collection	IgG EIA # 1
1 st	positive

This patient would be considered as already being immune.
4 weeks later:

sample collection	IgG EIA #1	IgG EIA #2	IgM ISAGA
1 st	positive	negative	negative
2 nd	positive	positive	positive

In fact, this patient is seroconverting. False positive or over-sensitivity of the EIA 1 test? Without the second sample collection, the EIA 1 test (IgG detected before IgM - rare) would have wrongly indicated this patient as being immunized, contrary to the EIA 2 test.

The true value of an IgG test is specificity for increased diagnostic accuracy.

B

Suspected infection ToRC IgM

Case study: Comparison of CMV IgM tests⁽⁷⁾

sample collection	IgM EIA #1, #3 and #4	IgM EIA #2	IgG
1 st	positive or equivocal	negative	positive

This patient would be suspected as being recently infected. In this example, the EIA 2 test appears to be less sensitive than the other tests (EIA 1, 3 and 4).

sample collection	IgG Avidity measurement
1 st	high avidity


High avidity excludes recent infection and therefore enables IgM to be interpreted as residual. In this example, the EIA 2 test avoided unnecessary patient distress and additional sample collections. The EIA 2 test detects less residual IgM.

The added value of an IgM test is specificity for less unnecessary distress.

C

Exclusion of recent infection IgG Avidity

Case study: Comparison of Toxo IgM tests with a Toxo Avidity test on a population initially found to be positive with the EIA 1 IgM test⁽⁸⁾



The true value of an Avidity test is exclusion of recent infection for less unnecessary investigations.

(a) 13 % and 24 % of the samples with the EIA 2 and EIA 3 tests respectively were not confirmed IgM positive (less residual IgM detection).

(b) 85 % of the samples presented a high avidity index, and recent infection was excluded.

VIDAS[®] ToRC PANEL



	VIDAS [®] ToRC	
VIDAS TOXO IgM	30202	60 tests
VIDAS TOXO IgG II	30210	60 tests
VIDAS TOXO IgG Avidity	30222	30 tests
VIDAS TOXO Competition	30211	60 tests
VIDAS CMV IgM	30205	30 tests
VIDAS CMV IgG	30204	60 tests
VIDAS CMV IgG Avidity II	413557	30 tests
VIDAS RUB IgM	30214	30 tests
VIDAS RUB IgG II	30221	60 tests

Other complementary tests

Toxo ISAGA	75361	192 tests
Toxo ISAGA IgA	79322	48 tests

REFERENCES

- Jenum PA, Stray-Pedersen B, Gundersen AG. **Improved diagnosis of primary *Toxoplasma gondii* infection in early pregnancy by determination of antitoxoplasma immunoglobulin G avidity.** *J Clin Microbiol.* 1997;35:1972-7.
- Wilson M, Remington JS, Clavet C, et al. **Evaluation of six commercial kits for detection of human immunoglobulin M antibodies to *Toxoplasma gondii*.** The FDA Toxoplasmosis Ad Hoc Working Group. *J Clin Microbiol.* 1997;35:3112-15.
- Grangeot-Keros L, Simon B, Audibert F, Vial M. **Should we routinely screen for cytomegalovirus antibody during pregnancy ?** *Intervirolgy* 1998;41:158-162.
- Denis F (Ed). **Les virus transmissibles de la mère à l'enfant.** John Libbey Eurotext 1999 (Pages 345-64).
- Ambroise-Thomas P, Petersen E (Eds.) **Congenital toxoplasmosis.** Springer 2000 (Pages 237-49).
- De Champs C, Pelloux H, Cambon M, et al. **Evaluation of the second generation IMx Toxo IgG antibody assay for detection of antibodies to *Toxoplasma gondii* in human sera.** *J Clin Lab Analysis* 1997;11:214-19.
- Fassotte R., Jost I. **Détection des IgM anti-cmv par technique Meia (AxSym, Abbot) à l'aide d'une nouvelle trousse de réactif à base d'antigène recombinant ; comparaison de trois autres méthodes de routine (VIDAS, Enzygnost, IMx).** *Immunoanal Biol Spéc.* 1999;14:260-63.
- van Gool T, Gilis H. **Avidity IgG : an effective diagnostic tool for unmasking false-positive IgM in toxoplasma serology.** Abstract presented during the European Conference on Congenital Toxoplasmosis, Vienna, June 29 - July 1, 2000.