





LUNDI 18 & MARDI 19 DÉCEMBRE 2023





P-166

Cytomegalovirus-specific T cell immunity: Design of an innovative functional assay

DANIEL Soizic¹, TRUFFOT Aurélie², CABEZAS Lara³, DUCROT Sandrine¹, RENARD Nathalie¹, PRETIS Catherine¹, MOREL Mathieu⁴, GERMI Raphaele², BERTHIER Franck⁵

¹R&D ImmunoAssays bioMérieux, Marcy l'Etoile, FRANCE; ² Univ. Grenoble Alpes, CNRS, CEA, IRIG IBS, Grenoble, FRANCE; ³ Univ. Grenoble Alpes, CNRS, Inserm, CHU Grenoble Alpes IAB, Grenoble, FRANCE; ⁴ Industrialisation Raw Materials bioMérieux, Marcy l'Étoile, FRANCE; ⁵ LifeSciences bioMérieux, Marcy l'Étoile, FRANCE

Franck.berthier@biomerieux.com - Tel: 06 26 88 73 86

Introduction

- Cytomegalovirus (CMV) responsible of 10-30 % of infections in kidney transplant recipients (Tedesco Silva Junior et al. Transpl Infect Dis. 2023)
- In immunocompromised patients: need to control immune system to avoid opportunistic infections and graft loss.
- Prevention: prophylactic treatment, follow up of CMV replication (RT-PCR), and evaluation of cell-mediated immunity (CMI) (Quantiferon, ELISpot) (Rutger Callens et al. Transplant Cell Ther. 2022)

Objective

- <u>Problematic</u>: Existing cellular diagnostic tools are very cumbersome and/or are limited to CD8+ T cell response monitoring.
- Need new assays (simpler, more automated and robust) to evaluate CMV-specific T cell and to stratify individual risk of CMV disease and then personalized prophylaxis

Evaluation of a new Interferon Gamma Release Assay (IGRA) on VIDAS ®.


A standardized and fully automated solution to facilitate quantification of CMI response to CMV on VIDAS®.

Method

1/ Workflow and analytical validation on fresh whole blood on healthy volunteers



2/ VIDAS® IGRA vs in house ELISpot T IFNy, on PBMC (at Grenoble University Hospital)



by VIDAS®, in supernatant

CMV specific lymphocyte T response, by in house

ELISpot T IFNy (Mabtech®)



* VIDAS® STIMM™ CMV RUO: Multiple CMV antigens at low concentration + immuno-booster designed to improved performances on immunocompromised patients

Results

1/ VIDAS® STIMM™ CMV RUO: Correlation between CMV T cell immunity and serology




Fig.1- VIDAS® STIMM™
CMV RUO presents a very good specificity and a good correlation with serology.

The two seroconverted patients were found positive with the IGRA test, highlighting the sensitivity of the technic.

3/ VIDAS® STIMM™ CMV RUO-IGRA: Towards a fully automated and standardized assay





Fig.3- 8 healthy donors CMV-IgG seropositive stimulated with VIDAS® STIMM™ CMV RUO optimized formulation in manual or fully automated version on VIDAS®3.

The fully automated version of the IGRA test contributes to high sensitivity in keeping good specificity

2/ VIDAS® STIMM™ CMV RUO: Impact of the formulation



- **Fig.2-** Whole blood of 3 CMV-IgG negative and 8 CMV-IgG positive stimulated with our CMV-specific antigens with or without optimized formulation.
- → Optimized formulation allows to increase specifically IFN-γ production by CMV-specific T cells

4/ In house ELISpot T IFNγ vs CMV-IGRA response on PBMC

- ➤ Positive ELISpot T result: >10 spots per well (250 000 cells)
- Positive VIDAS® CMV-IGRA result: IFNγ > 0,10 IU/mL and NIL < 0,08 IU/mL.</p>




Fig.4- Comparison of PBMC ELISpot T IFNy response with VIDAS® STIMM™ CMV RUO and others CMV (JPT) T-cell activated peptides.

Better correlation (R=0,84) were obtained with the JPT **pool pp65 + IE-1**

→ PP65 + IE-1 **Tab.1-** Comparison of qualitative results between VIDAS® STIMM™ CMV RUO (fluorescence value) vs ELISpot T IFNγ (RSV) Better concordance when using PP-65 or JPT pool (pp65 + IE-1)

Conclusion

- · A fully automated and standardized VIDAS® CMV-IGRA assay to detect CMV-specific cell mediated immunity in whole blood has been developed
- A sensitive and specific detection of cellular immune response to CMV in immunocompromised patients, with a unique CMV-specific antigen formulation BMX®
- Activation by VIDAS® STIMM™ CMV RUO correlate with stimulation by pp65 + IE-1 (JPT)
- Excellent concordance of CMV-IGRA assay (VIDAS® STIMM™ CMV RUO) and *in house* ELISpot with pp65 (JPT)