

BIOFIRE® SPOTFIRE® REDEFINING CONNECTIVITY



SEAMLESS CONNECTIVITY

To streamline your workflow, the SPOTFIRE System offers bidirectional Point-of-Care Testing (POCT) Connectivity to provide end-to-end data flow between your connected SPOTFIRE devices and your data management system.

ENHANCED CONNECTIVITY. SEAMLESS INTEROPERABILITY.

POCT Connectivity

Connect one or more SPOTFIRE® devices to your data manager via the industry-standard CLSI POCT01-A2 communication protocol to facilitate communication with your data manager in the background.

The SPOTFIRE System requires access to a local network and supports ethernet and Wi-Fi connectivity.

Quality Control Reporting

Automatically transmit Quality Control (QC) results to your data manager to help monitor and maintain test reliability across a network of BIOFIRE SPOTFIRE devices.

FASTER RESULTS. FEWER ERRORS.

Get Results Sooner

The SPOTFIRE System enables electronic reporting of patient test results, reducing manual entry and decreasing result turnaround time

Reduce Reporting Errors

The SPOTFIRE System accurately translates patient test results into a machine-readable XML format that eliminates the need for manual data entry and reduces the potential for misreported patient results.

MANAGEMENT ANYWHERE. ANYTIME.

Remotely Manage Operators

The SPOTFIRE System allows you to maintain control and oversight of operators and their training certifications at different locations across your organization.

Remotely Manage Device Settings

Configure QC requirements, electronic reporting options, operator login credentials, device date/time settings, and more.



CUSTOMER SUPPORT

Our goal is to provide you with the best support available while minimizing any disruption to your workflow. For BIOFIRE products technical assistance and support, please contact your local bioMérieux representative or distributor.

Contact Us

bioMérieux S.A.
69280 Marcy l'Etoile
France
Tel: +33 (0) 4 78 87 20 00
Fax: +33 (0) 4 78 87 20 90
biomerieux.com

Manufactured by:
BioFire Diagnostics, LLC
515 Colorow Drive
Salt Lake City, UT 84108 USA
Tel: +1-801-736-6354
biofiredx.com

VISIT BIOFIREDX.COM/COMPANY/CONTACT/
FOR MORE INFORMATION.

Product availability varies by country. Please contact your local bioMérieux representative for details.