



BIOMÉRIEUX

VIDAS[®] 3 for LIFE SCIENCES

Empower your research



PIONEERING DIAGNOSTICS

Meet any request at any time...

One of your priorities is to generate robust data.

VIDAS® brand is trusted around the world as a leader in on-demand immunoassay testings.

VIDAS® 3 empowers your research.

VIDAS® 3 is designed with your specific needs and workflow in mind to provide test results that are precise and accountable.

The system is equipped with automation, traceability and continual access, for improved lab efficiency.



Height: 61 cm; 24 in / Width: 75 cm; 29.5 in /
Depth: 64 cm; 25.5 in

VIDAS® 3



VIDAS® 3 allows fast adoption

Intuitive automation: test loading is streamlined, ensuring standardization within your lab

Fully automated pipettor action with single-use, disposable tips, to avoid contamination



VIDAS® 3 handles test security

Less human manipulation from the sample to the result

Sample and reagent barcodes are scanned upon test loading for complete traceability of all test components

Automated calibration, on-board dilution

Connected



VIDAS® 3 focuses on lab efficiency

Facilitated lab workflow with a simple answer to complex needs

Reduced operator time for maintenance, to better focus on high-level tasks

Remote access



VIDAS® 3 is designed for continual access

Three separate racks each containing three configurable segments can hold a total of **27 samples, diluents and calibrators** in the temperature-controlled loading bay

Four independent compartments can each run 3 tests



> VIDAS® 3 to perform your fully automated functional assays

- **A range of tests to explore cell-mediated immunity:** open IGRA assays, stimulants, cytokines
- **All steps to perform functional assays are done on board and fully automated:** pipeting, stimulation, incubation read-out.
- **Ready-to-use reagents**
- **One order reference per test:** kits contain all components needed for one test



VIDAS® 3 FOR LIFE SCIENCES



9325572 008/GBA This document is not legally binding. bioMérieux reserves the right to modify specifications without notice / BIOMÉRIEUX, the BIOMÉRIEUX logo, PIONEERING DIAGNOSTICS, VIDAS, are used, pending and/or registered trademarks belonging to bioMérieux or one of its subsidiaries or one of its companies. Any other name or trademark is the property of its respective owner / bioMérieux S.A. RCS Lyon 673 620 399 / Printed in France / EKNO / 537 519 449 RCS Lyon